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METTLER **TOLEDO**

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► www.mt.com/webinars

Petrochemicals
Analytical solutions in the laboratory

14
News
Edition 01/2012

Improved Water Analysis of Oils Combining Two Techniques

Critical mechanical systems, such as the wings of an aircraft, brakes in a car or pistons in an engine can corrode or malfunction if the lubricants and oil contain too much water. Quantitative water determination is a critical activity during the development and production of these substances. Various analytical methods have been developed over time to cope with this issue. The newest method is actually the past revisited.

Long before the ascent of Karl Fischer titration today's water determination method of choice. Edward Dean and David Stark developed a method using azeotropic distillation with toluene utilizing the fact that toluene is immiscible with water and also has a lower density and a higher boiling point than water.

Historic Dean-Stark distillation
Oil used to be weighed and mixed with toluene in a Dean-Stark apparatus. Subsequent heating caused the water to evaporate along with some toluene. Condensation and the subsequent separation of toluene and water happened in a graduated flask so that analysts could read the volume and calculate the water content of the original sample. The Dean-Stark method is still used today to determine substantial amounts of water (typically no less than 1%), but is inaccurate due to its reliance on visual water volume determination. It is also very time consuming, requiring a high level of manual interaction and tolerance.

Advantages of Karl Fischer
Many lubricants contain additives making them unsuitable for direct analysis in a Karl Fischer titration due to side reactions. Used engine oils can be directly analyzed but require frequent reagent and chemical changes wasting time and money. Drying ovens, in conjunction with a Karl Fischer titration, offer a solution but require high temperatures, are often time consuming and can experience tube blockages due to the evaporation of other volatiles, which can ultimately lead to system failure.

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GWPP[®] Recommendation: MTLABTEC-104-20090917-58

Balance
 IP2000N

Weighing Parameters
 Typical maximum weight for 0.1 % weighing accuracy (approximate factor $n = 3$, without safety factor)
 Tare: 20 mg Required: 50 mg
 Safety factor: Required: 3
 Allowed: 0.8 Required: 3
 For the determination of the actual minimum weight the issuing of a minimum weight certificate according to LSP is recommended.

Installation, qualification and training⁽¹⁾
 6 Month Plan

Service contract
 Calibration, Maintenance & Repair

Control of weighing equipment
 6000 Certificates⁽²⁾
 The routine testing functionalities of the IP Balance ensure that your balance is always tested correctly and on time.

Tests⁽¹⁾

Calibration (by service)	Quarterly
Repeatability (by user)	Monthly
Repeatability (by user)	Quarterly
Sensitivity (by user)	Monthly
FAT (by customer)	Yearly

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For correct handling of weights we recommend the use of appropriate tongs and/or gloves.
 Recalibration interval of weights: Weekly test plan

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© 04/2012 Mettler-Toledo AG